


Title: <i>In planta</i> co-immunoprecipitation (Summer School protocol)	
Reference No: TSL037	Version No: 1.0

THE SAINSBURY LABORATORY STANDARD OPERATING PROCEDURE		<small>TheSainsburyLaboratory</small> 
TITLE: <i>In planta</i> co-immunoprecipitation (Summer School protocol)		
APPLIES TO STAFF IN: Sainsbury Laboratory and Students at the summer school		
HEALTH & SAFETY INFORMATION INCLUDED: YES		
REFERENCE No: TSL037	VERSION No: 1.0	
DATE EFFECTIVE: August 2017	REVIEW DATE: August 2019	
AUTHOR: Joe Win	APPROVED BY: Simon Foster	

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: <i>In planta</i> co-immunoprecipitation (Summer School protocol)	
Reference No: TSL037	Version No: 1.0

1 PURPOSE OF PROCEDURE/METHOD AND ITS SCOPE

To purify plant proteins that associate with the effectors expressed *in planta*. *Nicotiana benthamiana* will be used as a model host plant.

2 EQUIPMENT NEEDED

Pipettes
Pipette tips
Centrifuge
Anti-FLAG M2 affinity gel
GTEN: 10% glycerol, 25 mM Tris pH 7.5, 1 mM EDTA, 500 mM NaCl
Extraction Buffer: GTEN, 2% w/v PVPP, 10 mM DTT, 1x protease inhibitor cocktail (Sigma), 0.5% Igepal
Immunoprecipitation (IP) buffer: GTEN, 0.5% Igepal
Conservation Buffer: 5 mM Tris pH 7.5
Frozen leaf powder obtained from *N. benthamiana* leaves expressing the effectors
Dry ice

3 STEPS IN PROCEDURE

- **Wear lab coat, safety glasses, and gloves. Observe local safety rules.**
- **Be careful not to make long contact with dry ice when handling the tube with frozen leaf powder**
- **Dispose of all biologically contaminated waste according to local safety rules, e. g., use the 'kill box' for re-useable items and blue box for disposable items. Sharp items must be enclosed in yellow-cap plastic bottles.**

- 3.1 Weigh out 1 g leaf powder and resuspend it in 2 ml ice-cold extraction buffer. Vortex thoroughly until the solution is homogenous, and keep on ice.
- 3.2 Centrifuge at full speed in a microcentrifuge for 10 min at 4°C. Transfer the supernatant into syringe and filter through 0.45 µm membranes. Keep the filtrate (~1 ml) in a new tube and place on ice.
- 3.3 Re-suspend the anti-FLAG gel by tapping the side of the vial several times and mix with a cut pipette tip.

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: <i>In planta</i> co-immunoprecipitation (Summer School protocol)	
Reference No: TSL037	Version No: 1.0

- 3.4 Pipette enough gel (20 μ L per sample) into a 1.5 ml Eppendorf tube. Centrifuge at $800 \times g$ for 1 min and carefully remove the supernatant (**take care not to aspirate the gel**).
- 3.5 Re-suspend the gel in 1 ml of IP buffer.
- 3.6 Centrifuge at $800 \times g$ for 1 min and remove the supernatant.
- 3.7 Re-suspend the gel to 2x original volume with IP buffer and add 40 μ L of diluted gel to the leaf extract prepared above (step 2).
- 3.8 Incubate the gel and the leaf extract at 4°C for 1-3 h on a rotating mixer.
- 3.9 Centrifuge at $800 \times g$ for 30 s. Discard supernatant and add 1 ml of fresh IP buffer. Repeat four more times but always leave about 50 μ l at the bottom of the tube to avoid aspirating the beads. After the last wash, remove as much supernatant as possible without touching the beads. Beads are now ready for trypsin digestion and mass spectrometry.




4 RISK STATEMENT

This activity is low risk and does not involve handling any hazardous chemicals except EDTA at very low levels (1 mM) and small quantities of DTT and IGEPAL CA-630.

All individuals using this procedure will be shown the risk assessment and given appropriate information, instruction and training in the risks and precautions necessary, including the use of any personal protective equipment required.

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

SOP HEALTH RISK ASSESSMENT

[1] Activity:	<i>In planta</i> co-immunoprecipitation				
[2] Location of activity:	The Sainsbury Laboratory and Chris Lamb Training Suite				
[3] Who is involved:	Science staff in The Sainsbury Laboratory and participants to the TSL Summer School				
[4] Frequency of activity:	Variable				
[5] Duration of Activity:	7-10 days				
[6] Chemical Hazard Name:	Hazard Statements	Route of Exposure*	Quantity Used		
Ethylenediaminetetraacetic acid (EDTA) 	H319- Causes serious eye irritation	S, I	Not more than 2 ml at 1 mM		
Dithiothreitol (DTT) 	H302- Harmful if swallowed H315- Causes skin irritation H319- Causes serious eye irritation H412- Harmful to aquatic life with long lasting effects	S, I	<1g		
IGEPAL® CA-630 	H302- Harmful if swallowed H315- Causes skin irritation H318- Causes serious eye damage H410- Very toxic to aquatic life with long lasting effects	S, I	<1 ml		
[7] Details of biological agents of risk to human health:			GMRA Number		
None					
[8] Other Hazards: Please <input checked="" type="checkbox"/> as necessary					
Hot or Cold Burns	<input checked="" type="checkbox"/>	Ionising Radiation**	<input type="checkbox"/>	Ultra Violet or Infra Red	<input type="checkbox"/>
Dust	<input type="checkbox"/>	Noise	<input type="checkbox"/>	Pollen Sensitizer	<input type="checkbox"/>

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Repetitive Action	<input type="checkbox"/>	Extreme Cold Environment (< 0°C)	<input type="checkbox"/>	Lifting / Manual Handling	<input type="checkbox"/>
Asphyxiation	<input type="checkbox"/>	Cuts	<input type="checkbox"/>	Electrical	<input type="checkbox"/>
Slips / trips / falls	<input type="checkbox"/>	Display Screen Equipment	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>
[9] Control Measures: Please <input type="checkbox"/> as necessary					
Fume Cupboard	<input type="checkbox"/>	Microbiological Safety Cabinet	<input type="checkbox"/>	Total Containment Cabinet	<input type="checkbox"/>
Ventilated Bench	<input type="checkbox"/>	Spill Tray	<input type="checkbox"/>	Trained personnel only	<input checked="" type="checkbox"/>
Signs	<input type="checkbox"/>	Reduce frequency/alternate activity	<input type="checkbox"/>	Reduce duration of activity	<input type="checkbox"/>
Sub divide a load	<input type="checkbox"/>	2 man lift of equipment	<input type="checkbox"/>	Not for more than 1 hour	<input type="checkbox"/>
Regular, short breaks	<input type="checkbox"/>	Alternate activities	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>
[10] Personal Protection: Please <input type="checkbox"/> as necessary					
Lab coat	<input checked="" type="checkbox"/>	Safety Glasses	<input checked="" type="checkbox"/>	Face Shield	<input type="checkbox"/>
Goggles	<input type="checkbox"/>	Gloves	<input checked="" type="checkbox"/>	Thermal Protective Gloves	<input type="checkbox"/>
Ear defenders	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>		
[11] Is personal monitoring and/or health surveillance required? Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>					
Details:					
[12] Restrictions: Please <input type="checkbox"/> as necessary					
No lone working	<input type="checkbox"/>	Not to be left unattended	<input type="checkbox"/>	Named persons only	<input type="checkbox"/>
In restricted area	<input type="checkbox"/>	Not by new or expectant mothers	<input type="checkbox"/>	Under constant supervision	<input type="checkbox"/>
Not by under 18's	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>		
[13] Level of Residual Risk: Please <input type="checkbox"/> as necessary					
Low	<input checked="" type="checkbox"/>	Medium	<input type="checkbox"/>	High	<input type="checkbox"/>

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Name of Assessor:

Date:

Simon Foster

21st July 2017

* Route of exposure; S = skin, I = ingestion, B = inhalation

Activities involving new or expectant mothers and young persons also require additional risk assessment.

1 DOCUMENTATION

Links to relevant H&S information on intranet or internet

Reference any relevant manuals

Link to: JIC Chemical Tables:

http://intranet/infoserv/support/QualityAssurance/Chemical_Tables_SOPs.htm

Link to: Good Laboratory Practice in the Use of Chemicals:

http://intranet/infoserv/support/Safety/Chemical/GLP_Chems.htm

Link to Biological and GM Safety:

<http://intranet/infoserv/Support/Safety/Biological/index.htm>

Link to Laboratory Waste Disposal:

<http://intranet/infoserv/Support/Safety/Waste/index.htm>

Win J et al., (2011) Purification of Effector–Target Protein Complexes via Transient Expression in *Nicotiana benthamiana*. In *Plant Immunity*, McDowell JM (ed) pp 181-194-194. Humana Press

Petre B et al., (2015) Candidate Effector Proteins of the Rust Pathogen *Melampsora larici-populina* Target Diverse Plant Cell Compartments. *Mol Plant Microbe Interact* 28: 689-700

6 RELATED PROCEDURES

Other relevant SOPs

7 NOTES

8 APPENDICES

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: <i>In planta</i> co-immunoprecipitation (Summer School protocol)	
Reference No: TSL037	Version No: 1.0

Chemical	Ethylenediaminetetraacetic acid (EDTA)	Sigma: 431788	CAS: 60-00-4
Hazard Statement	H319 Causes serious eye irritation.		
Precautionary Statement	P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.		
Handling	Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.		
Storage	Store in cool place. Keep container tightly closed in a dry and well-ventilated place.		
Disposal	Dispose of in a chemically compatible container; ensure liquids are placed in a container designed for liquids. Label and place in the collection tray for "Waste Not Suitable for Bulking-Up" in the Chemical Waste Store. Where possible, use the container in which the chemical was supplied. Follow the user instructions for Using the Chemical Waste Store . Consult the Chemical Safety Officer or your Lab Manager for more information.		
Spillage	Wear appropriate personal safety equipment; lab coat, gloves and safety glasses and absorb spill with inert material e.g. vermiculite and place in suitable container and dispose of in Chemical Waste Store.		

Chemical	Dithiothreitol (DTT)	Sigma Aldrich D0632	CAS:3483-12-3
Hazard Statement(s)	H302 Harmful if swallowed. H315 Causes skin irritation. H319 Causes serious eye irritation.		
Precautionary Statement(s)	P301 + P312 + P330 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.		
Handling	Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed.		
Storage	Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Recommended storage temperature: 2 - 8 °C		
Disposal	Dispose of in a chemically compatible container; ensure liquids are placed in a container designed for liquids. Label and place in the collection tray for "Waste Not Suitable for Bulking-Up" in the Chemical Waste Store. Where possible, use the container in which the chemical was supplied. Follow the user instructions for Using the Chemical Waste Store . Consult the Chemical Safety Officer or your Lab Manager for more information.		
Spillage	Wear appropriate personal safety equipment; lab coat, gloves and safety glasses and clean up spills with blue roll or spill tamer kit. Put in suitable container and dispose of in Chemical Waste Store.		

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used.
This copy was printed on 21 July 2017.

Title: <i>In planta</i> co-immunoprecipitation (Summer School protocol)	
Reference No: TSL037	Version No: 1.0

Chemical	Igepal® CA-630	Sigma I8896	CAS: 9002-93-1
Hazard Statement(s)	H302 Harmful if swallowed. H315 Causes skin irritation. H318 Causes serious eye damage. H410 Very toxic to aquatic life with long lasting effects. .		
Precautionary Statement(s)	P280 Wear eye protection/ face protection. P301 + P312 + P330 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.		
Handling	Avoid contact with skin and eyes. Avoid inhalation of vapour or mist		
Storage	Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.		
Disposal	Dispose of in a chemically compatible container; ensure liquids are placed in a container designed for liquids. Label and place in the collection tray for "Waste Not Suitable for Bulking-Up" in the Chemical Waste Store. Where possible, use the container in which the chemical was supplied. Follow the user instructions for Using the Chemical Waste Store . Consult the Chemical Safety Officer or your Lab Manager for more information.		
Spillage	Wear appropriate personal safety equipment; lab coat, gloves and safety glasses and clean up spills with blue roll or spill tamer kit. Put in suitable container and dispose of in Chemical Waste Store.		

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used.
This copy was printed on 21 July 2017.