


Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

THE SAINSBURY LABORATORY STANDARD OPERATING PROCEDURE		<small>TheSainsburyLaboratory</small> 
TITLE: Yeast-two-Hybrid assay (Summer School Protocol)		
APPLIES TO STAFF IN: Sainsbury Laboratory		
HEALTH & SAFETY INFORMATION INCLUDED: YES		
REFERENCE No: TSL039	VERSION No: 1.0	
DATE EFFECTIVE: August 2017	REVIEW DATE: August 2019	
AUTHOR: Joe Win	APPROVED BY: Simon Foster	

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

1 PURPOSE OF PROCEDURE/METHOD AND ITS SCOPE

To test protein-protein interaction in yeast.

2 EQUIPMENT NEEDED

Pipettes
 Pipette tips
 Water bath
 Shaking incubator
 Y2H Gold cells from Clontech
 Cuvettes
 1M lithium acetate (LiAc) solution
 50% PEG (MW. 3350) solution
 1M TE-buffer pH 8
 0.2% adenine solution
 Dimethyl sulfoxide (DMSO)
 Carrier-DNA (for transformation)
 Plasmid DNA
 Yeast Dropout Medium (-leu/-trp and -leu/-trp/-ade/-his)
 YPD-Medium

3 STEPS IN PROCEDURE

This yeast two hybrid (Y2H) protocol is based on the user manual of Clontech MatchMaker GAL4 Two Hybrid System 3. Please consult this manual for further details and trouble-shooting.

- **Wear lab coat and gloves. Observe local safety rules. Do not eat or drink in the laboratory.**
- **Dispose of all biologically contaminated waste according to local safety rules, e. g., use the 'kill box' for re-useable items and blue box for disposable items. Sharp items must be enclosed in yellow-cap plastic bottles.**

Yeast Transformation

1. Grow Y2H Gold cells in 50 ml YPD (+Adenine) liquid medium in a shaking incubator at 28 °C at 220 RPM overnight up to OD₆₀₀ 1.0
2. Spin cells at 1100 x g for 5 min in a centrifuge at room temperature

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

3. Wash the cells with H₂O twice
4. Resuspend the cells in 5 ml sterile 100 mM LiAc, 100mM TE-buffer pH 8.0 solution
5. Mix 100 ng each of two plasmid DNAs (in 1-5 µl volume) with 10 µl of carrier DNA
6. Add 100 µl cells to the DNA mix
7. Add 500 µl sterile LiAc/PEG solution (100 mM LiAc; 40% PEG, 100 mM TE-buffer pH 8.0), mix well and incubate in a shaking heat block at 28 °C for 30 min
8. Add 70 µl sterile DMSO, mix gently
9. Incubate at 42 °C for 15 min in a water bath, cool on ice for 2 min
10. Spin cells at 16000 x g) for 30 sec at room temperature
11. Resuspend cells in 250 µl 100 mM TE-buffer, plate on SD (-leu/-trp) agar plates

Y2H assay

1. Inoculate 3 ml SD -leu/-trp liquid culture with a single transformant (single colony) from the SD -leu/-trp selection plate, and place at 28 °C shaking incubator at 220 RPM overnight
2. Measure the OD₆₀₀ of the overnight culture (dilute culture 1:10 in SD medium for measurement)
3. Adjust the OD₆₀₀ of all test strains to 1.0
4. Prepare serial 1:10 dilutions (1:10 / 1:100 / 1:1000 i.e. OD₆₀₀ 1 / 0.1 / 0.01 / 0.001). Mix well by pipetting
5. Prepare replica drop plates of all dilutions on SD -leu/-trp and SD -leu/-trp/-ade/-his+α-X-gal by placing 3 µl drops from each dilution. Let the droplets dry on the plate
6. Incubate the plates at in an incubator at 28 °C for 4 days
7. Record the results by taking a photograph/scan of each plate

4 RISK STATEMENT


This activity is low risk and does not involve handling of any hazardous chemicals except adenine which is used in growing media at trace amount.

All individuals using this procedure will be shown the risk assessment and given appropriate information, instruction and training in the risks and

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

precautions necessary, including the use of any personal protective equipment required.

SOP HEALTH RISK ASSESSMENT				
[1] Activity:	Yeast-two-hybrid assay to test protein-protein interaction			
[2] Location of activity:	The Sainsbury Laboratory and Chris Lamb Training Suite			
[3] Who is involved:	Science staff in The Sainsbury Laboratory and participants to the TSL Summer School			
[4] Frequency of activity:	Variable			
[5] Duration of Activity:	3 hours			
[6] Chemical Hazard Name:	Hazard Statements	Route of Exposure*	Quantity Used	
Lithium acetate	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.	N/A	Less than 1 ml	
Polyethylene glycol	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.	N/A	Less than 1 ml	
Dimethyl sulfoxide	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.	N/A	Less than 1 ml	
Adenine 	H301- Toxic if swallowed.	I	Less than 1 ml	
[7] Details of biological agents of risk to human health:			GMRA Number	
None				
[8] Other Hazards: Please <input checked="" type="checkbox"/> as necessary				
Hot or Cold Burns	<input type="checkbox"/>	Ionising Radiation**	<input type="checkbox"/>	Ultra Violet or Infra Red <input type="checkbox"/>

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

Dust	<input type="checkbox"/>	Noise	<input type="checkbox"/>	Pollen Sensitizer	<input type="checkbox"/>		
Repetitive Action	<input type="checkbox"/>	Extreme Cold Environment (< 0°C)	<input type="checkbox"/>	Lifting / Manual Handling	<input type="checkbox"/>		
Asphyxiation	<input type="checkbox"/>	Cuts	<input type="checkbox"/>	Electrical	<input type="checkbox"/>		
Slips / trips / falls	<input type="checkbox"/>	Display Screen Equipment	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>		
[9] Control Measures: Please <input type="checkbox"/> as necessary							
Fume Cupboard	<input type="checkbox"/>	Microbiological Safety Cabinet	<input type="checkbox"/>	Total Containment Cabinet	<input type="checkbox"/>		
Ventilated Bench	<input type="checkbox"/>	Spill Tray	<input type="checkbox"/>	Trained personnel only	<input checked="" type="checkbox"/>		
Signs	<input type="checkbox"/>	Reduce frequency/alternate activity	<input type="checkbox"/>	Reduce duration of activity	<input type="checkbox"/>		
Sub divide a load	<input type="checkbox"/>	2 man lift of equipment	<input type="checkbox"/>	Not for more than 1 hour	<input type="checkbox"/>		
Regular, short breaks	<input type="checkbox"/>	Alternate activities	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>		
[10] Personal Protection: Please <input type="checkbox"/> as necessary							
Lab coat	<input checked="" type="checkbox"/>	Safety Glasses	<input type="checkbox"/>	Face Shield	<input type="checkbox"/>		
Goggles	<input type="checkbox"/>	Gloves	<input checked="" type="checkbox"/>	Thermal Protective Gloves	<input type="checkbox"/>		
Ear defenders	<input type="checkbox"/>	Other (give details)	<input type="checkbox"/>				
[11] Is personal monitoring and/or health surveillance required? Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>							
Details:							
[12] Restrictions: Please <input type="checkbox"/> as necessary							
No lone working	<input type="checkbox"/>	Not to be left unattended	<input type="checkbox"/>	Named persons only	<input type="checkbox"/>		
In restricted area	<input type="checkbox"/>	Not by new or expectant mothers	<input type="checkbox"/>	Under constant supervision	<input type="checkbox"/>		
Not by under 18's	<input checked="" type="checkbox"/>	Other (give details)	<input type="checkbox"/>				
[13] Level of Residual Risk: Please <input type="checkbox"/> as necessary							
Low	<input checked="" type="checkbox"/>	Medium	<input type="checkbox"/>	High	<input type="checkbox"/>		

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

Name of Assessor:	Date:
Simon Foster	21 st July 2017

* Route of exposure; S = skin, I = ingestion, B = inhalation

Activities involving new or expectant mothers and young persons also require additional risk assessment.

5 DOCUMENTATION

Links to relevant H&S information on intranet or internet

Reference any relevant manuals

Link to: JIC Chemical Tables:

http://intranet/infoserv/support/QualityAssurance/Chemical_Tables_SOPs.htm

Link to: Good Laboratory Practice in the Use of Chemicals:

http://intranet/infoserv/support/Safety/Chemical/GLP_Chems.htm

Link to Biological and GM Safety:

<http://intranet/infoserv/Support/Safety/Biological/index.htm>

Link to Laboratory Waste Disposal:

<http://intranet/infoserv/Support/Safety/Waste/index.htm>

6 RELATED PROCEDURES

Other relevant SOPs

7 NOTES

8 APPENDICES

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used. This copy was printed on 21 July 2017.

Title: Yeast-two-Hybrid assay	
Reference No: TSL039	Version No: 1.0

Chemical	Lithium acetate dihydrate	Sigma L6883	CAS: 6108-17-4
Hazard Statement(s)	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008. This substance is not classified as dangerous according to Directive 67/548/EEC.		
Precautionary Statement(s)	N/A		
Handling	Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection.		
Storage	Store in cool place. Keep container tightly closed in a dry and well-ventilated place.		
Disposal	Dispose of in a chemically compatible container; ensure liquids are placed in a container designed for liquids. Label and place in the collection tray for "Waste Not Suitable for Bulking-Up" in the Chemical Waste Store. Where possible, use the container in which the chemical was supplied. Follow the user instructions for Using the Chemical Waste Store . Consult the Chemical Safety Officer or your Lab Manager for more information.		
Spillage	Wear appropriate personal safety equipment; lab coat, gloves and safety glasses and clean up spills with blue roll or spill tamer kit. Put in suitable container and dispose of in Chemical Waste Store.		

Chemical	Dimethyl Sulphoxide (DMSO)	Sigma D2650	CAS:67-68-5
Hazard Statement(s)	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.		
Precautionary Statement(s)	N/A.		
Handling	Avoid inhalation of vapour or mist. Keep away from sources of ignition - No smoking. Take measures to prevent the build up of electrostatic charge.		
Storage	Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Store under inert gas. Hygroscopic		
Disposal	Dispose of in a chemically compatible container; ensure liquids are placed in a container designed for liquids. Label and place in the collection tray for "Waste Not Suitable for Bulking-Up" in the Chemical Waste Store. Where possible, use the container in which the chemical was supplied. Follow the user instructions for Using the Chemical Waste Store . Consult the Chemical Safety Officer or your Lab Manager for more information.		
Spillage	Wear appropriate personal safety equipment; lab coat, gloves and safety glasses and clean up spills with blue roll or spill tamer kit. Put in suitable container and dispose of in Chemical Waste Store.		

This is a controlled document maintained on electronic media. When appearing in paper form it should be checked against the master on the QA database to ensure that the current version is being used.
This copy was printed on 21 July 2017.