



Yeast two hybrid: Basics

- Method to investigate protein-protein interaction
 - Based on the transcription factor Gal4
 - Gal4 consists of a DNA-binding domain (BD) and an activation domain (AD)
 - Both domains are functional even when separated
 - Both domains need to be in close proximity to activate transcription
-
- Generate fusion proteins in autonomously replicating plasmids (1)
 - Fusion proteins: 1st POI-Gal4-BD (Bait) / 2nd POI-Gal4-AD (Prey)
 - Both plasmids are co-transformed into Gal4-deficient Yeast strain (Y2H Gold; (2))
 - Plasmids complement Leucine and Tryptophan auxotrophy (3)
(Selection on SD –leu/-trp medium)
-
- Bait fusion protein can bind to upstream activating sequences (UAS) of Gal4-responsive promoter regions
 - Interacting proteins bring BD and AD in close proximity
 - POI1-BD/POI2-AD complex activate transcription of reporter genes (4)
(Histidine / Adinine / α -Galactosidase / Aureobasidine)
(selection on SD –leu/-trp/-his/-ade + α -X-Gal)



Yeast two hybrid: The two main methods (Library Screen)

Library screen



Your Gene of interest

Sequence AD-Plasmid
(to identify interaction partner)






Figure 5. Screening an AD fusion library using strain AH109. Use the stringency of your choice to screen for interacting proteins. Note: high stringency selections result in fewer colonies, and reduce the number of false positives. However, weak interactions may be missed.

From: Matchmaker™ GAL4 Two-Hybrid System 3 & Libraries User Manual

Yeast two hybrid: The two main methods (Co-transformation)



Material provided:

- Y2H Gold o/n culture (in YPDA medium)
 - Solutions to prepare competent cells (H₂O, Te/LiAc)
 - Solutions for Yeast transformation (PEG/LiAc/TE, DMSO)
 - Plasmids for Yeast Transformation (3x bait plasmid containing target genes, 2x prey plasmid containing effectors)
-
- Test strains o/n culture (in SD 2XDO medium)
 - Control strains o/n culture (in SD 2XDO medium)
 - Y2H Media (2x DO, 4x DO Medium)

Yeast two hybrid: Pros and Cons

Pros:

- Quick and easy screening method
- Yeast is easy to handle in the lab (growth, transformation, selection)
- Semi *in-vivo* conditions
- Allows for post-translational modifications (e.g. Glycosylation, Prenylation etc.)
- Yeast chaperones support proper protein folding

Cons:

- Depends on nuclear localization
 - Post-translational modifications in yeast can be different to the natural system
 - Co-factors might not be present in yeast
 - False negatives:
 - Wrong fold of proteins can alter interaction
 - TM-domains or Prenylation can interfere with nuclear localization
 - False positives:
 - Wrong fold can cause unnatural interaction
 - Proteins which are naturally separated (spatial or temporal) can interact in yeast
 - Intrinsic activation domain can cause autoactivation in BD-fusion proteins
- Needs downstream verification by additional experiments

Yeast two hybrid: Verification of Interactions

Methods to verify interaction partners:

- FRET (Förster resonance energy transfer)
- Co-Immunoprecipitation (CoIP)
- *in-vitro* interaction studies (e.g. analytical gel filtration)

Pros:

- *in planta* conditions, transient expression, does not require protein extraction
- *in planta* conditions, transient expression, quick
- “clean” environment (no protein modifications/additional proteins), detection of binding affinity, allows to study the mode of interaction (e.g. by protein structure resolution)

Additional Yeast based assays:

- | | |
|--------------------------|---|
| - Yeast-one-Hybrid | - screen for Protein-DNA interaction |
| - Yeast-three-hybrid | - screen for Protein-RNA interaction |
| - Split-ubiquitin system | - screen for Protein-Protein interaction of membrane proteins |